


Product Datasheet

Malondialdehyde antibody [Out of stock] (orb396641)


Description

Mouse monoclonal antibody against


Species/Host	Mouse
Reactivity	All
Conjugation	ATTO-390
Tested Applications	ELISA, FACS, FC, ICC, WB
Immunogen	Synthetic Malondialdehyde modified Keyhole Limpet Hemocyanin (KLH).
Target	Malondialdehyde
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% Sodium azide *Storage buffer may change when conjugated
Form/Appearance	PBS pH 7.4, 50% glycerol, 0.09% Sodium azide *Storage buffer may change when conjugated
Concentration	p>1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:1000 dilution of SMC-514 was sufficient for detection of Malondialdehyde in 2 µg of Malondialdehyde conjugated to BSA by ECL immunoblot analysis using Goat Anti-Mouse IgG:HRP as the secondary Antibody.
Isotype	IgG1
Clonality	Monoclonal
Purity	Protein G Purified
Dilution Range	WB (1:1000); ICC/IF (1:50); FACS (1:50); FCM (1:50); ELISA (1:1000)
Expiration Date	12 months from date of receipt.



Flow Cytometry analysis of human Neurobl...



Western Blot analysis of Lane 1: Molecula...



Immunofluorescence analysis of human Emb...