



Product Datasheet

Malondialdehyde Antibody (orb396658)


Description	Mouse Anti-Malondialdehyde (MDA) Monoclonal IgG1
Species/Host	Mouse
Reactivity	All
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IHC, WB
Immunogen	Synthetic Malondialdehyde modified Keyhole Limpet Hemocyanin (KLH).
Target	Malondialdehyde
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% Sodium azide *Storage buffer may change when conjugated
Form/Appearance	PBS pH 7.4, 50% glycerol, 0.09% Sodium azide
Concentration	1 mg/ml
Storage	-20°C
Note	For research use only
Application notes	A 1:1000 dilution of SMC-515 was sufficient for detection of Malondialdehyde in 2 µg of Malondialdehyde conjugated to BSA by ECL immunoblot analysis using Goat Anti-Mouse IgG:HRP as the secondary Antibody.
Isotype	IgG1
Clonality	Monoclonal
Purity	Protein G Purified
Dilution Range	WB (1:1000); ICC/IF (1:50); ELISA (1:1000); IHC-P (1:100)
Expiration Date	12 months from date of receipt.



Western Blot analysis of Human Cervical ...



Western Blot analysis of Lane 1: Molecul...



Immunocytochemistry analysis of human Em...