


Product Datasheet

Rat TROY ELISA Kit (orb397145)

Description	Sandwich ELISA for quantitative detection of Rat TROY concentrations.
Reactivity	Rat
Range	62.5 pg/ml - 4000 pg/ml
Tested Applications	ELISA
Target	TNFRSF19
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Detection Range: Cell Culture Supernates, Serum, Plasma, Cell Lysates, Tissue Homogenates
Sample Types	Cell Culture Supernates, Serum, Plasma, Cell Lysates, Tissue Homogenates
Uniprot ID	D3ZUJ0
Sensitivity	30 pg/ml
Entrez	290300
Expiration Date	Please enquire.



Representative
Standard
Curve