

## Product Datasheet

### CD303 antibody (FITC) (orb397146)

**Description**

CD303 antibody (FITC)

**Species/Host**

Rabbit

**Reactivity**

Human, Mouse, Rat

**Conjugation**

FITC

**Tested**

ICC, IF

**Applications**
**Immunogen**

KLH conjugated synthetic peptide derived from human CD303.

**Target**

CD303 (FITC)

**Form/Appearance**

Liquid: PBS, 2 mM Sodium Azide, 30% glycerol

**Concentration**

- 100 µg (in 200 µl): 0.5 mg/ml

**Storage**

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

**Note**

For research use only

**Isotype**

IgG

**Clonality**

Polyclonal

**Purity**

Polyclonal antibodies are purified by peptide affinity chromatography

**MW**

25 kDa

**Uniprot ID**
[Q8WTT0](#)
**NCBI**
[NM\\_130441.2](#), [NM\\_130441.21](#),  
[NP\\_569708.1](#), [NP\\_569708.11](#)
**Entrez**
[170482](#)
**Dilution Range**

IF/ICC: 1:100-300

**Expiration Date**

12 months from date of receipt.