

Product Datasheet

CD163SR1 antibody (orb397158)

Description	CD163SR1 antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Protein
Target	CD163SR1
Form/Appearance	10 mM PBS, 0.02% sodium azide
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Affinity purified by Protein A
Expiration Date	12 months from date of receipt.