

Product Datasheet

CCL5 antibody (orb397164)

Description	CCL5 antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Protein
Target	CCL5
Form/Appearance	10 mM PBS, 0.02% sodium azide
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Affinity purified by Protein A
Expiration Date	12 months from date of receipt.