

Product Datasheet

SOS1 antibody (orb39718)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	SOS1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	152464
Uniprot ID	Q07889
NCBI	NP_005624.2
Dilution Range	WB 1:50-1:100
Expiration Date	12 months from date of receipt.

