



Product Datasheet

DLK2 antibody (orb39719)


Description	Rabbit polyclonal antibody to DLK2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	DLK2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	40548
Uniprot ID	Q6UY11
NCBI	NP_996262.1 , NP_076421.2
Dilution Range	WB: 1:1000, IHC-P: 1:50-100, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Western blot analysis of DLK2 antibody i...



Flow cytometric analysis of DLK2 antibody...