

Product Datasheet

Phospho-Smad3(S425) Mouse Monoclonal Antibody(1D9) (orb397216)

Description	Phospho-Smad3(S425) Mouse Monoclonal Antibody(1D9)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Monoclonal
Uniprot ID	P84022
Expiration Date	12 months from date of receipt.