

Product Datasheet

HGD antibody (orb39755)

Description	HGD antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Homogentisate 1,2-dioxygenase protein (226-445AA)
Target	HGD
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q93099
Dilution Range	IHC-P:1:20-200
Expiration Date	12 months from date of receipt.