

Product Datasheet

LEFTY2 antibody (orb39761)

Description	LEFTY2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Left-right determination factor 2 protein (57-366AA)
Target	LEFTY2
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O00292
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.