




Product Datasheet

UTY antibody (orb39763)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Synthetic Peptide
Target	UTY
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	149548
Uniprot ID	O14607
NCBI	NP_009056.3 , NP_001245183.1 , NP_872600.1 , NP_872601.1 , NP_001245181.1
Dilution Range	WB: 1:1000, FACS: 1:10-50, IF/ICC: 1:10-50
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of MDA-MB231...



Western blot analysis of in mouse brain ...