




Product Datasheet

PKM2 antibody (orb39765)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	PKM
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	57937
Uniprot ID	P14618
NCBI	NP_001193728.1 , NP_001193727.1 , NP_001193726.1 , NP_872270.1 , NP_872271.1 , NP_002645.3 , NP_001193725.1
Dilution Range	WB: 1:1000, IHC-P: 1:50-100, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis
of 293 cell lysate...



Immunohistochemical
staining of paraffin...