




Product Datasheet

EFEMP1 antibody (orb39786)

Description	Rabbit polyclonal antibody to EFEMP1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	EFEMP1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	54641
Uniprot ID	Q12805
NCBI	NP_001034437.1 , NP_001034438.1
Dilution Range	WB: 1:1000, IHC-P: 1:50-100, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Flow cytometry analysis of EFEMP1 antibody...



Immunohistochemical analysis of formalin...