



Product Datasheet

Aggrecan Antibody (orb398162)

Description	Aggrecan Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Concentration	1.0mg/mL
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P16112
Dilution Range	WB: 1:500-1000, IF/ICC: 1:50-1:200, IHC-P: 1:50-1:100, ELISA: 1:10000-1:20000
Expiration Date	12 months from date of receipt.



Western blot analysis of Aggrecan expres...



Immunohistochemical analysis of Aggrecan...