

Product Datasheet

GABARAPL1 antibody (orb39824)

Description Peptide affinity purified rabbit polyclonal antibody.

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications FC, IHC-P, WB

Immunogen Recombinant Protein

Target GABARAPL1

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal



MW 14044

Uniprot ID [Q9H0R8](#)


NCBI [NP_113600.1](#)

Dilution Range WB: 1:1000, IHC-P: 1:50-100, FACS: 1:10-50

Expiration Date 12 months from date of receipt.



Immunohistochemical staining of paraffin...



Western blot analysis of in HeLa cell li...