

Product Datasheet

5-hydroxymethylcytosine antibody (orb39852)

Description	Rabbit polyclonal antibody to 5-hydroxymethylcytosine is vital in genome function. It is dynamically...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Target	This 5-hydroxymethylcytosine antibody is generated from rabbits immunized with human 5-hydroxymethylcytosine.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.