



Product Datasheet

KIF3A antibody (orb39854)


Description	KIF3A antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Kinesin-like protein KIF3A protein (430-699AA)
Target	KIF3A
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9Y496
Dilution Range	WB:1:1000-5000, IHC-P:1:20-200
Expiration Date	12 months from date of receipt.



Western blot aanalysis
of MCF-7 whole ce...



Immunohistochemistry
of paraffin-embedde...



Immunohistochemistry
of paraffin-embedde...