


Product Datasheet

MBTPS1 antibody (orb39861)

Description	MBTPS1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human Membrane-bound transcription factor site-1 protease protein (803-1052AA)
Target	MBTPS1
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q14703
Dilution Range	IHC-P:1:20-200
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse kidney tissue showing a single band at approximately 85kDa, indicating the presence of MBTPS1 protein.