

Product Datasheet

FLAD1 antibody (orb39892)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	FLAD1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	65266
Uniprot ID	Q8NFF5
NCBI	NP_079483.3 , NP_001171821.1 , NP_001171820.1 , NP_958800.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

