




Product Datasheet

TSPAN6 antibody (orb39919)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthetic Peptide
Target	TSPAN6
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	27563
Uniprot ID	O43657
NCBI	NP_001265672.1 , NP_003261.1
Dilution Range	WB 1:50-1:100
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
HeLa cells u...



Western Blot
of m.bladder
cell lysate
wi...