

Product Datasheet

Aquaporin 1 antibody (orb399652)

Description	Aquaporin 1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	FITC
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human Aquaporin 1.
Target	Aquaporin 1
Form/Appearance	Liquid: PBS, 2 mM Sodium Azide, 30% glycerol
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	29 kDa
Uniprot ID	Q02013, P29975, P29972
Entrez	358
Dilution Range	IF/ICC: 1:100-300
Expiration Date	12 months from date of receipt.