

Product Datasheet

BAZ2A antibody (FITC) (orb399665)

Description	BAZ2A antibody (FITC)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	FITC
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human BAZ2A.
Target	BAZ2A (FITC)
Form/Appearance	Liquid: PBS, 2 mM Sodium Azide, 30% glycerol
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	211 kDa
Uniprot ID	Q9UIF9 , Q91YE5
NCBI	NP_038477.2 , NM_013449.3
Entrez	11176
Dilution Range	IF/ICC: 1:100-300
Expiration Date	12 months from date of receipt.