

Product Datasheet

Beclin 1 antibody (FITC) (orb399667)

Description

Beclin 1 antibody (FITC)

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	FITC
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human Beclin 1.
Target	Beclin 1 (FITC)
Form/Appearance	Liquid: PBS, 2 mM Sodium Azide, 30% glycerol
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	52 kDa
Uniprot ID	O88597 , Q91XJ1 , Q14457
NCBI	NP_003757.1 , NM_003766.31 , NP_003757.11 , NM_003766.3
Entrez	8678
Dilution Range	IF/ICC: 1:100-300
Expiration Date	12 months from date of receipt.