

Product Datasheet

CD38 antibody (FITC) (orb399702)

Description

CD38 antibody (FITC)

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

FITC

Tested

ICC, IF

Applications
Immunogen

KLH conjugated synthetic peptide derived from human CD38.

Target

CD38 (FITC)

Form/Appearance

Liquid: PBS, 2 mM Sodium Azide, 30% glycerol

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

34 kDa

Uniprot ID
P56528, Q64244
NCBI
NP_031672.21, NP_031672.2, NM_007646.41, NM_007646.4
Entrez
12494
Dilution Range

IF/ICC: 1:100-300

Expiration Date

12 months from date of receipt.