

Product Datasheet

COL2A1 antibody (FITC) (orb399723)

Description

COL2A1 antibody (FITC)

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

FITC

Tested

ICC, IF

Applications
Immunogen

KLH conjugated synthetic peptide derived from human COL2A1.

Target

COL2A1 (FITC)

Form/Appearance

Liquid: PBS, 2 mM Sodium Azide, 30% glycerol

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

117 kDa

Uniprot ID
P02458, P05539, P28481
NCBI
NP_001835.31, NM_001844.4, NM_001844.41, NP_001835.3
Entrez
1280
Dilution Range

IF/ICC: 1:100-300

Expiration Date

12 months from date of receipt.