

Product Datasheet

Estrogen Receptor beta 2 antibody (FITC) (orb399764)

Description

Estrogen Receptor beta 2 antibody

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	FITC
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human Estrogen Receptor beta 2.
Target	Estrogen Receptor beta 2 (FITC)
Form/Appearance	Liquid: PBS, 2 mM Sodium Azide, 30% glycerol
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	59 kDa
Uniprot ID	P11475 , Q61539 , O95718
NCBI	NM_004452.3 , NP_004443.3
Entrez	2103
Dilution Range	IF/ICC: 1:100-300
Expiration Date	12 months from date of receipt.