

Product Datasheet

Gremlin 1 antibody (FITC) (orb399783)

Description	Gremlin 1 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	FITC
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human Gremlin 1.
Target	Gremlin 1 (FITC)
Form/Appearance	Liquid: PBS, 2 mM Sodium Azide, 30% glycerol
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	21 kDa
Uniprot ID	070326, O60565, O35793
Entrez	26585
Dilution Range	IF/ICC: 1:100-300
Expiration Date	12 months from date of receipt.