

## Product Datasheet

### MAP2K7 antibody (FITC) (orb399828)

|                            |                                                                                                                |
|----------------------------|----------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | MAP2K7 antibody (FITC)                                                                                         |
| <b>Species/Host</b>        | Rabbit                                                                                                         |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                              |
| <b>Conjugation</b>         | FITC                                                                                                           |
| <b>Tested Applications</b> | ICC, IF                                                                                                        |
| <b>Immunogen</b>           | KLH conjugated synthetic peptide derived from human MAP2K7.                                                    |
| <b>Target</b>              | MAP2K7 (FITC)                                                                                                  |
| <b>Form/Appearance</b>     | Liquid: PBS, 2 mM Sodium Azide, 30% glycerol                                                                   |
| <b>Concentration</b>       | - 100 µg (in 200 µl): 0.5 mg/ml                                                                                |
| <b>Storage</b>             | Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles. |
| <b>Note</b>                | For research use only                                                                                          |
| <b>Isotype</b>             | IgG                                                                                                            |
| <b>Clonality</b>           | Polyclonal                                                                                                     |
| <b>Purity</b>              | Polyclonal antibodies are purified by peptide affinity chromatography                                          |
| <b>MW</b>                  | 47 kDa                                                                                                         |
| <b>Dilution Range</b>      | IF/ICC: 1:100-300                                                                                              |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                |