

Product Datasheet

SLC22A8 antibody (FITC) (orb399922)

Description	SLC22A8 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	FITC
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human SLC22A8.
Target	SLC22A8 (FITC)
Form/Appearance	Liquid: PBS, 2 mM Sodium Azide, 30% glycerol
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	60 kDa
Uniprot ID	Q8TCC7
Entrez	9376
Dilution Range	IF/ICC: 1:100-300
Expiration Date	12 months from date of receipt.