

Product Datasheet

Syndecan 2 antibody (FITC) (orb399936)

Description	Syndecan 2 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	FITC
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human Syndecan 2.
Target	Syndecan 2 (FITC)
Form/Appearance	Liquid: PBS, 2 mM Sodium Azide, 30% glycerol
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	22 kDa
Uniprot ID	P34741, P34900, P43407
Entrez	6383
Dilution Range	IF/ICC: 1:100-300
Expiration Date	12 months from date of receipt.