



Product Datasheet

Sgta antibody (orb400103)

Description	Goat polyclonal antibody to Sgta
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	Sgta (mouse aa84-97)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 35kDa band observed in lysates of cell line PD19 and doublet of 30+35kDa in lysates of cell line K562 (calculated MW of 34.3kDa according to NP_078775.1). Recommended concentration: 1-2
Protein Sequence	PSEEDSAEAERLKT
Clonality	Polyclonal
MW	34.3
Entrez	6449
Dilution Range	ELISA: 1:2000, WB: 1-2 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Testis tissue showing Sgta expression.



Western blot analysis of PD19 (lane 1) and K562 (lane 2) cell lysates showing Sgta protein bands.