



Product Datasheet

Lmo3 antibody (orb400111)

Description	Goat polyclonal antibody to Lmo3
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC
Target	LIM domain only 3 / Lmo3
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Immunofluorescence: Customer results showed positive staining in Lmo3heterozygous mutant mice (Publication adding on this page very soon...).
Protein Sequence	EGLMKEGYAPQVR
Clonality	Polyclonal
MW	16.6; 17.9; 18.4; 19.0
Entrez	55885
Dilution Range	ELISA: 1:128000, IHC-P: 10ug/ml
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of paraforma...



Flow cytometric analysis of paraformalde...