

Product Datasheet

CDK9 antibody (orb400114)

Description	Goat polyclonal antibody to CDK9
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	CDK9
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 40-45kDa band observed in nuclear cell lysates of NIH3T3 (calculated MW of 42.8kDa according to Mouse NP_570930.1). Recommended concentration: 0.3-1
Protein Sequence	KGSQITQQSTNQSRN
Clonality	Polyclonal
MW	42.8
Entrez	1025
Dilution Range	ELISA: 1:8000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.