

Product Datasheet

mpt51 antibody (FITC) (orb400158)

Description	mpt51 antibody (FITC)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	FITC
Immunogen	Recombinant Mycobacterium tuberculosis MPT51/MPB51 antigen protein (27-299AA)
Target	mpt51
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P9WQN7
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.