


Product Datasheet

DLX6 antibody (orb400192)

Description	DLX6 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Homeobox protein DLX-6 protein (110-171AA)
Target	DLX6
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P56179
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293 whole ce...