

Product Datasheet

I79_001071 antibody (FITC) (orb400266)

Description	I79_001071 antibody (FITC)
Species/Host	Rabbit
Reactivity	Hamster
Conjugation	FITC
Immunogen	Recombinant Cricetulus griseus C-C motif chemokine protein (25-143AA)
Target	I79_001071
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	G3GTT2
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.