



Product Datasheet

HLA-DRB3 antibody (orb400272)

Description	HLA-DRB3 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human HLA class II histocompatibility antigen, DR beta 3 chain protein (111-227AA)
Target	HLA-DRB3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P79483
Dilution Range	WB:1:500-5000, IF:1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of Raji whole cell...



Immunofluorescence analysis of A549 cell...