



Product Datasheet

COPS3 antibody (orb40031)


Description	COPS3 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Immunogen	Recombinant Human COP9 signalosome complex subunit 3 protein (194-423AA)
Target	COPS3
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UNS2
Dilution Range	WB:1:1000-5000, IHC-P:1:20-200, IP:1:200-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa(lane 1), 2...



Immunohistochemical staining of human te...



Immunoprecipitation analysis of HeLa who...