

Product Datasheet

ARPC5L antibody (Biotin) (orb400347)

Description	ARPC5L antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Actin-related protein 2/3 complex subunit 5-like protein (15-124AA)
Target	ARPC5L
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9BPX5
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.