



Product Datasheet

GTPBP4 antibody (orb400380)

Description	GTPBP4 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human Nucleolar GTP-binding protein 1 protein (347-634AA)
Target	GTPBP4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9BZE4
Dilution Range	WB:1:500-5000, IF:1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293 whole ce...



Immunofluorescence analysis of HeLa cell...