



Product Datasheet

ACAA2 antibody (orb400392)


Description	ACAA2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	Recombinant Human 3-ketoacyl-CoA thiolase, mitochondrial protein (188-257AA)
Target	ACAA2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P42765
Dilution Range	WB:1:1000-5000, IHC-P:1:500-1000, IF:1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of Hela whole cell...



Immunohistochemical staining of Human ki...



Immunofluorescence analysis of HepG2 cel...