



Product Datasheet

CAMLG antibody (orb400416)


Description	CAMLG antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IP, WB
Immunogen	Recombinant Human Calcium signal-modulating cyclophilin ligand protein (1-189AA)
Target	CAMLG
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P49069
Dilution Range	WB:1:500-5000, IF:1:50-500, IP:1:200-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat whole ce...



Immunofluorescence analysis of HepG2 cel...



Immunoprecipitation analysis of HEK293 W...