



Product Datasheet

CCNH antibody (orb400420)

Description	CCNH antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Cyclin-H protein (261-319AA)
Target	CCNH
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P51946
Dilution Range	IHC-P:1:20-200
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse liver tis...



Immunohistochemical staining of Human ad...