



Product Datasheet

GCNT3 antibody (orb400456)


Description	GCNT3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Beta-1,3-galactosyl-O-glycosyl-glycoprotein beta-1,6-N-acetylglucosaminyltransferase 3 protein (250-419AA)
Target	GCNT3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O95395
Dilution Range	IHC-P:1:20-200
Expiration Date	12 months from date of receipt.



Western blot analysis of A549 whole cell...



Immunohistochemical staining of human pa...



Immunohistochemical staining of human co...