


Product Datasheet

KIF20A antibody (orb400488)

Description	KIF20A antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP
Immunogen	Recombinant Human Kinesin-like protein KIF20A protein (161-294AA)
Target	KIF20A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O95235
Dilution Range	IHC-P:1:50-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human te...