

Product Datasheet

RPS27A antibody (FITC) (orb400522)

Description	RPS27A antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Ubiquitin-40S ribosomal protein S27a protein (16-56AA)
Target	RPS27A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P62979
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.