

Product Datasheet

TBL1X antibody (Biotin) (orb400559)

Description	TBL1X antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human F-box-like/WD repeat-containing protein TBL1X protein (1-300AA)
Target	TBL1X
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O60907
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.