


Product Datasheet

yebF antibody (orb400612)

Description	yebF antibody
Species/Host	Rabbit
Reactivity	E. coli
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Escherichia coli Protein YebF protein (22-118AA)
Target	yebF
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P33219
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.



Western blot analysis of Recombinant pro...