


Product Datasheet

TEF antibody (orb400628)

Description	TEF antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Thyrotroph embryonic factor protein (7-303AA)
Target	TEF
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q10587
Dilution Range	WB:1:500-5000
Expiration Date	12 months from date of receipt.



Western blot analysis of U87 whole cell ...